



SLIDE WARMING TABLE

SW-HEAT-PRO-24

Design, Development & Manufacturing of Histopathology Laboratory Instruments.



9034999314



ipisales63@gmail.com



ipihisto.com



Proudly made in



SLIDE WARMING TABLE



SW-HEAT-PRO-24

Slide Warming Table is used in histopathology laboratory for drying slides at a uniform temperature just below the melting point of wax so as to provide flattening without damage to the section. slide warming table is fabricated from top notch materials as per defined industry standards. Our Slide Warming Table are known for their optimum utility and high performance.

KEY FEATURES



CONSTANT TEMPERATURE

Maintains uniform and stable temperature for consistent performance.



TEFLON COATED TOP

Ensures smooth surface, corrosion resistance and easy cleaning.



SILICON HEATER

Provides efficient and uniform heating for reliable operation.



USER FRIENDLY

Simple to operate with intuitive controls for convenience.



MICROPROCESSOR CONTROL

Advanced control for precise temperature regulation and enhanced accuracy.



DIGITAL DISPLAY

Clear and accurate readout for easy monitoring.

TECHNICAL SPECIFICATION:

SPECIFICATION	DETAILS
Display	LCD display showing: Date, Time, Day, Set value, Present value.
Temperature Control	Microprocessor Precise Temperature control with accuracy.
Uniform Temperature	Uniform temperature throughout the surface.
Construction	Resistant material for durability.
Temperature Range	Ambient to 100°C
Temperature Tolerance	+/- 1°C
Slide Warmer Dimension	240 x 240 mm
Power	220V 50Hz
Overall Dimension	360 x 260 x 125 mm



9034999314



ipisales63@gmail.com




ipihisto.com



Proudly made in





SLIDE WARMING TABLE

SW-HEAT-PRO-24

Design, Development & Manufacturing of Histopathology Laboratory Instruments.

~ INTERNATIONAL PATHALOGY INSTRUMENTS ~



9034999314



ipisales63@gmail.com



ipihisto.com



Proudly made in

