



PARAFFIN WAX DISPENSER

WAX-PARA-24

Design, Development & Manufacturing of Histopathology Laboratory Instruments.



9034999314



ipisales63@gmail.com



ipihisto.com



Proudly made in





PARAFFIN WAX DISPENSER

WAX-PARA-24

Paraffin Wax Bath (**WAX-PARA-24**) is used to melt solid paraffin wax, maintain it in liquid form, and dispense it whenever required. It prevents rust formation and features a temperature-controlled nozzle to ensure smooth wax flow and avoid blockage caused by solid or semi-solid wax.



KEY FEATURES

- Stainless steel body with rust-resistant design
- Fast and uniform wax melting system
- Thermostatic temperature control
- Maintains wax in molten form for continuous use
- Temperature-controlled nozzle for smooth wax flow
- Prevents nozzle blockage due to solidified wax
- Easy dispensing and drainage system
- Safe, durable, and user-friendly operation

TECHNICAL SPECIFICATION

- **Body Material:** Stainless Steel
- **Heating Type:** Electric Heating Element
- **Temperature Range:** 40°C to 80°C (adjustable)
- **Temperature Control:** Thermostatic Controller
- **Power Supply:** 220-240V AC, 50/60 Hz
- **Power Consumption:** 1-2 kW
- **Dispensing System:** Temperature-controlled nozzle/tap
- **Capacity** - 2, 5, 10 litre

